

TNFRSF10D Antibody, HRP conjugated



PACO32557

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for the cytotoxic ligand TRAIL. Contains a truncated death domain and hence is not capable of inducing apoptosis but protects against TRAIL-mediated apoptosis. Reports are contradictory with regards to its ability to induce the NF-kappa-B pathway. According to Ref.1, it cannot but according to Ref.2, it can induce the NF-kappa-B pathway.

Gene ID:

TNFRSF10D

Uniprot

Q9UBN6

Synonyms:

Tumor necrosis factor receptor superfamily member 10D (Decoy receptor 2) (DcR2) (TNF-related apoptosis-inducing ligand receptor 4) (TRAIL receptor 4) (TRAIL-R4) (TRAIL receptor with a truncated death domain) (CD antigen CD264), TNFRSF10D, DCR2 TRAILR4 TRUNDD

Immunogen:

Recombinant Human Tumor necrosis factor receptor superfamily member 10D protein (58-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A