

Ccl4 Antibody



PACO32656

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Monokine with inflammatory and chemokinetic properties.

Gene ID:

Ccl4

Uniprot

P14097

Synonyms:

C-C motif chemokine 4 (Immune activation protein 2) (ACT-2) (ACT2) (Macrophage inflammatory protein 1-beta) (MIP-1-beta) (Protein H400) (SIS-gamma) (Small-inducible cytokine A4), Ccl4, Mip1b Scya4


Immunogen:

Recombinant Mouse C-C motif chemokine 4 protein (24-92AA).

Storage:

Preservative: 0.1mol/L NaCl, 15mmol/L NaN₃, PH7.2

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO32656 at dilution of 1:100.