

PACO32660

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Monokine with inflammatory, pyrogenic and chemokinetic properties. Circulates at high concentrations in the blood of healthy animals. Binding to a high-affinity receptor activates calcium release in neutrophils. It also inhibits colony formation of bone marrow myeloid immature progenitors.

Gene ID:

Ccl9

Uniprot

P51670

Synonyms:

C-C motif chemokine 9 (CCF18) (Macrophage inflammatory protein 1-gamma) (MIP-1-gamma) (Macrophage inflammatory protein-related protein 2) (MRP-2) (Small-inducible cytokine A9) [Cleaved into: CCL9(29-101); CCL9(30-101); CCL9(31-101)], Ccl9, Mrp2 Scya10 Scya9

Immunogen:

Recombinant Mouse C-C motif chemokine 9 protein (22-114AA).

Storage:

Preservative: 0.1mol/L NaCl, 15mmol/L NaN3, PH7.2

Product Images

N/A

N/A