

PACO32676

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Ligand for CXCR2. Has chemotactic activity for neutrophils. May play a role in inflammation and exert its effects on endothelial cells in an autocrine fashion.

Gene ID:

Cxcl3

Uniprot

Q6W5C0

Synonyms:

C-X-C motif chemokine 3 (Dendritic cell inflammatory protein 1), Cxcl3, Dcip1 Gm1960

Immunogen:

Recombinant Mouse C-X-C motif chemokine 3 protein (32-100AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A