
Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be a cytokine that affects the growth, movement, or activation state of cells that participate in immune and inflammatory response.

Gene ID:

Cxcl9

Uniprot

P18340

Synonyms:

C-X-C motif chemokine 9 (Gamma-interferon-induced monokine) (Monokine induced by interferon-gamma) (MIG) (MuMIG) (Protein m119) (Small-inducible cytokine B9), Cxcl9, Mig Scyb9

Immunogen:

Recombinant Mouse C-X-C motif chemokine 9 protein (22-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A