

TACR1 Antibody, HRP conjugated



PACO32721

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is a receptor for the tachykinin neuropeptide substance P. It is probably associated with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinity of this receptor to tachykinins is: substance P > substance K > neuromedin-K.

Gene ID:

TACR1

Uniprot

P25103

Synonyms:

Substance-P receptor (SPR) (NK-1 receptor) (NK-1R) (Tachykinin receptor 1), TACR1, NK1R TAC1R

Immunogen:

Recombinant Human Substance-P receptor protein (1-29AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A