

PACO32732

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Flagellin is the subunit protein which polymerizes to form the filaments of bacterial flagella.

Gene ID:

fliC

Uniprot

P04949

Synonyms:

Flagellin, fliC, flaF hag

Immunogen:

Recombinant Escherichia coli Flagellin protein (2-498AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A