

# ISG15 Antibody, Biotin conjugated



PACO32739

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ISG15

**Uniprot**

P05161

**Synonyms:**

Ubiquitin-like protein ISG15 (Interferon-induced 15 kDa protein) (Interferon-induced 17 kDa protein) (IP17) (Ubiquitin cross-reactive protein) (hUCRP), ISG15, G1P2 UCRP

**Immunogen:**

Recombinant Human Ubiquitin-like protein ISG15 protein (2-157AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A