

WFDC2 Antibody, Biotin conjugated



PACO32759

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

WFDC2

Uniprot

Q14508

Synonyms:

WAP four-disulfide core domain protein 2 (Epididymal secretory protein E4) (Major epididymis-specific protein E4) (Putative protease inhibitor WAP5), WFDC2, HE4 WAP5

Immunogen:

Recombinant Human WAP four-disulfide core domain protein 2 protein (31-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A