

GDF15 Antibody, Biotin conjugated



PACO32807

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Growth differentiation factor-15 (GDF-15) is disclosed along with its polynucleotide sequence and amino acid, sequence. Also disclosed are diagnostic and therapeutic methods of using the GDF-15 polypeptide and polynucleotide sequences.

Gene ID:

GDF15

Uniprot

Q99988

Synonyms:

Growth/differentiation factor 15 (GDF-15) (Macrophage inhibitory cytokine 1) (MIC-1) (NSAID-activated gene 1 protein) (NAG-1) (NSAID-regulated gene 1 protein) (NRG-1) (Placental TGF-beta) (Placental bone morphogenetic protein) (Prostate differentiation factor), GDF15, MIC1 PDF PLAB PTGFB

Immunogen:

Recombinant Human Growth/differentiation factor 15 protein (198-308AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A