

# PYGB Antibody, HRP conjugated



PACO32817

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Glycogen phosphorylase that regulates glycogen mobilization. Phosphorylase is an important allosteric enzyme in carbohydrate metabolism. Enzymes from different sources differ in their regulatory mechanisms and in their natural substrates. However, all known phosphorylases share catalytic and structural properties.

**Gene ID:**

PYGB

**Uniprot**

P11216

**Synonyms:**

Glycogen phosphorylase, brain form (EC 2.4.1.1), PYGB

**Immunogen:**

Recombinant Human Glycogen phosphorylase, brain form protein (6-843AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A