

# GOLM1 Antibody, Biotin conjugated



PACO32879

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Unknown. Cellular response protein to viral infection.

**Gene ID:**

GOLM1

**Uniprot**

Q8NBJ4

**Synonyms:**

Golgi membrane protein 1 (Golgi membrane protein GP73) (Golgi phosphoprotein 2), GOLM1, C9orf155 GOLPH2

**Immunogen:**

Recombinant Human Golgi membrane protein 1 protein (36-401AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A