

# MUC16 Antibody



PACO32880

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:10-1:100

**Protein Background:**

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

**Gene ID:**

MUC16

**Uniprot**

Q8WXI7

**Synonyms:**

Mucin-16 (MUC-16) (Ovarian cancer-related tumor marker CA125) (CA-125) (Ovarian carcinoma antigen CA125), MUC16, CA125

**Immunogen:**


Recombinant Human Mucin-16 protein (21712-21973AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human prostate tissue using PACO32880 at dilution of 1:20.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO32880 at dilution of 1:20.