

CTSS Antibody, FITC conjugated



PACO32902

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Thiol protease. Key protease responsible for the removal of the invariant chain from MHC class II molecules. The bond-specificity of this proteinase is in part similar to the specificities of cathepsin L.

Gene ID:

CTSS

Uniprot

P25774

Synonyms:

Cathepsin S (EC 3.4.22.27), CTSS

Immunogen:

Recombinant Human Cathepsin S protein (115-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A