

# LRG1 Antibody, HRP conjugated



PACO32953

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LRG1

**Uniprot**

P02750

**Synonyms:**

Leucine-rich alpha-2-glycoprotein (LRG), LRG1, LRG

**Immunogen:**

Recombinant Human Leucine-rich alpha-2-glycoprotein protein (36-346AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A