

AQP4 Antibody, Biotin conjugated



PACO32963

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Forms a water-specific channel. Osmoreceptor which regulates body water balance and mediates water flow within the central nervous system.

Gene ID:

AQP4

Uniprot

P55087

Synonyms:

Aquaporin-4 (AQP-4) (Mercurial-insensitive water channel) (MIWC) (WCH4), AQP4

Immunogen:

Recombinant Human Aquaporin-4 protein (253-323AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A