

# HAGH Antibody, Biotin conjugated



PACO32971

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Thiolesterase that catalyzes the hydrolysis of S-D-lactoyl-glutathione to form glutathione and D-lactic acid, .

**Gene ID:**

HAGH

**Uniprot**

Q16775

**Synonyms:**

Hydroxyacylglutathione hydrolase, mitochondrial (EC 3.1.2.6) (Glyoxalase II) (Glx II), HAGH, GLO2 HAGH1

**Immunogen:**

Recombinant Human Hydroxyacylglutathione hydrolase, mitochondrial protein (50-308AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A