

# KRT8 Antibody, FITC conjugated



PACO32986

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

**Gene ID:**

KRT8

**Uniprot**

P05787

**Synonyms:**

Keratin, type II cytoskeletal 8 (Cytokeratin-8) (CK-8) (Keratin-8) (K8) (Type-II keratin Kb8), KRT8, CYK8

**Immunogen:**

Recombinant Human Keratin, type II cytoskeletal 8 protein (2-483AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A