

PACO32992

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.

Gene ID:

PGA3

Uniprot

P0DJJ8

Synonyms:

Pepsin A-3 (EC 3.4.23.1) (Pepsinogen-3), PGA3

Immunogen:

Recombinant Human Pepsin A-3 protein (16-388AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A