

PGR Antibody, FITC conjugated



PACO33018

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

The steroid hormones and their receptors are involved in the regulation of eukaryotic gene expression and affect cellular proliferation and differentiation in target tissues. Progesterone receptor isoform B (PRB) is involved activation of c-SRC/MAPK signaling on hormone stimulation.

Gene ID:

PGR

Uniprot

P06401

Synonyms:

Progesterone receptor (PR) (Nuclear receptor subfamily 3 group C member 3), PGR, NR3C3

Immunogen:

Recombinant Human Progesterone receptor protein (4-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A