

COL4A1 Antibody, HRP conjugated



PACO33041

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

COL4A1

Uniprot

P02462

Synonyms:

Collagen alpha-1(IV) chain [Cleaved into: Arresten], COL4A1

Immunogen:

Recombinant Human Collagen alpha-1(IV) chain protein (30-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A