

# COL4A1 Antibody, FITC conjugated



PACO33042

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COL4A1

**Uniprot**

P02462

**Synonyms:**

Collagen alpha-1(IV) chain [Cleaved into: Arresten], COL4A1

**Immunogen:**

Recombinant Human Collagen alpha-1(IV) chain protein (30-167AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A