

PAX2 Antibody, HRP conjugated



PACO33069

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable transcription factor that may have a role in kidney cell differentiation. Has a critical role in the development of the urogenital tract, the eyes, and the CNS.

Gene ID:

PAX2

Uniprot

Q02962

Synonyms:

Paired box protein Pax-2, PAX2

Immunogen:

Recombinant Human Paired box protein Pax-2 protein (1-363AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A