

# PAX2 Antibody, FITC conjugated



PACO33070

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable transcription factor that may have a role in kidney cell differentiation. Has a critical role in the development of the urogenital tract, the eyes, and the CNS.

**Gene ID:**

PAX2

**Uniprot**

Q02962

**Synonyms:**

Paired box protein Pax-2, PAX2

**Immunogen:**

Recombinant Human Paired box protein Pax-2 protein (1-363AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A