

TGM2 Antibody, FITC conjugated



PACO33126

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

Gene ID:

TGM2

Uniprot

P21980

Synonyms:

Protein-glutamine gamma-glutamyltransferase 2 (EC 2.3.2.13) (Tissue transglutaminase) (Transglutaminase C) (TG(C)) (TGC) (TGase C) (Transglutaminase H) (TGase H) (Transglutaminase-2) (TGase-2), TGM2

Immunogen:

Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 protein (4-544AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A