

BHMT2 Antibody, Biotin conjugated



PACO33135

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the regulation of homocysteine metabolism. Converts homocysteine to methionine using S-methylmethionine (SMM) as a methyl donor.

Gene ID:

BHMT2

Uniprot

Q9H2M3

Synonyms:

S-methylmethionine--homocysteine S-methyltransferase BHMT2 (SMM-hcy methyltransferase) (EC 2.1.1.10) (Betaine--homocysteine S-methyltransferase 2), BHMT2

Immunogen:

Recombinant Human S-methylmethionine--homocysteine S-methyltransferase BHMT2 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A