

BGLAP Antibody, Biotin conjugated



PACO33139

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Constitutes 1-2% of the total bone protein. It binds strongly to apatite and calcium.

Gene ID:

BGLAP

Uniprot

P02818

Synonyms:

Osteocalcin (Bone Gla protein) (BGP) (Gamma-carboxyglutamic acid, containing protein), BGLAP

Immunogen:

Recombinant Human Osteocalcin protein (52-100AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A