

# TTR Antibody, FITC conjugated



PACO33166

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

**Gene ID:**

TTR

**Uniprot**

P02766

**Synonyms:**

Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB

**Immunogen:**

Recombinant Human Transthyretin protein (21-147AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A