

TTR Antibody, Biotin conjugated



PACO33167

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Gene ID:

TTR

Uniprot

P02766

Synonyms:

Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB

Immunogen:

Recombinant Human Transthyretin protein (21-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A