

INS Antibody, Biotin conjugated



PACO33171

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Insulin decreases blood glucose concentration. It increases cell permeability to monosaccharides, amino acid, and fatty acid, . It accelerates glycolysis, the pentose phosphate cycle, and glycogen synthesis in liver.

Gene ID:

INS

Uniprot

P01308

Synonyms:

Insulin [Cleaved into: Insulin B chain; Insulin A chain], INS

Immunogen:

Recombinant Human Insulin protein (25-110AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A