

APOC3 Antibody, Biotin conjugated



PACO33195

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

APOC3

Uniprot

P02656

Synonyms:

Apolipoprotein C-III (Apo-CIII) (ApoC-III) (Apolipoprotein C3), APOC3

Immunogen:

Recombinant Human Apolipoprotein C-III protein (21-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A