

TG Antibody, Biotin conjugated



PACO33219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Gene ID:

TG

Uniprot

P01266

Synonyms:

Thyroglobulin (Tg), TG

Immunogen:

Recombinant Human Thyroglobulin protein (2500-2767AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A