

PACO33220

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Stimulates development of follicle and spermatogenesis in the reproductive organs.

Gene ID:

FSHB

Uniprot

P01225

Synonyms:

Follitropin subunit beta (Follicle-stimulating hormone beta subunit) (FSH-B) (FSH-beta) (Follitropin beta chain), FSHB

Immunogen:

Recombinant Human Follitropin subunit β protein (19-129AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A