

MMP12 Antibody, Biotin conjugated



PACO33231

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in tissue injury and remodeling. Has significant elastolytic activity. Can accept large and small amino acid, at the P1' site, but has a preference for leucine. Aromatic or hydrophobic residues are preferred at the P1 site, with small hydrophobic residues (preferably alanine) occupying P3.

Gene ID:

MMP12

Uniprot

P39900

Synonyms:

Macrophage metalloelastase (MME) (EC 3.4.24.65) (Macrophage elastase) (ME) (hME) (Matrix metalloproteinase-12) (MMP-12), MMP12, HME

Immunogen:

Recombinant Human Macrophage metalloelastase protein (109-470AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A