

RNASE3 Antibody, FITC conjugated



PACO33258

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cytotoxin and helminthotoxin with low-efficiency ribonuclease activity. Possesses a wide variety of biological activities. Exhibits antibacterial activity, including cytoplasmic membrane depolarization of preferentially Gram-negative, but also Gram-positive strains. Promotes E.coli outer membrane detachment, alteration of the overall cell shape and partial loss of cell content.

Gene ID:

RNASE3

Uniprot

P12724

Synonyms:

Eosinophil cationic protein (ECP) (EC 3.1.27. -) (Ribonuclease 3) (RNase 3), RNASE3, ECP RNS3

Immunogen:

Recombinant Human Eosinophil cationic protein (28-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A