

ENO3 Antibody, HRP conjugated



PACO33305

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Appears to have a function in striated muscle development and regeneration.

Gene ID:

ENO3

Uniprot

P13929

Synonyms:

Beta-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (Enolase 3) (Muscle-specific enolase) (MSE) (Skeletal muscle enolase), ENO3

Immunogen:

Recombinant Human Beta-enolase protein (7-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A