

PAGE1 Antibody



PACO33316

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cellular defense response.

Gene ID:

PAGE1

Uniprot

O75459

Synonyms:

P antigen family member 1 (PAGE-1) (AL5) (G antigen 9) (GAGE-9) (G antigen family B member 1) (Prostate-associated gene 1 protein), PAGE1, GAGE9 GAGEB1


Immunogen:

Recombinant Human P antigen family member 1 protein (1-146AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO33316 at dilution of 1:100.