

PIAS2 Antibody, FITC conjugated



PACO33330

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PIAS2

Uniprot

O75928

Synonyms:

E3 SUMO-protein ligase PIAS2 (EC 2.3.2. -) (Androgen receptor-interacting protein 3) (ARIP3) (DAB2-interacting protein) (DIP) (E3 SUMO-protein transferase PIAS2) (Msx-interacting zinc finger protein) (Miz1) (PIAS-NY protein) (Protein inhibitor of activated STAT x) (Protein inhibitor of activated STAT2), PIAS2, PIASX

Immunogen:

Recombinant Human E3 SUMO-protein ligase PIAS2 protein (1-572AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A