

STEAP1 Antibody, HRP conjugated



PACO33333

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

STEAP1

Uniprot

Q9UHE8

Synonyms:

Metalloreductase STEAP1 (EC 1.16.1. -) (Six-transmembrane epithelial antigen of prostate 1), STEAP1, PRSS24 STEAP

Immunogen:

Recombinant Human Metalloreductase STEAP1 protein (3-69AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A