

# KLK10 Antibody, HRP conjugated



PACO33405

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Has a tumor-suppressor role for NES1 in breast and prostate cancer.

**Gene ID:**

KLK10

**Uniprot**

O43240

**Synonyms:**

Kallikrein-10 (EC 3.4.21. -) (Normal epithelial cell-specific 1) (Protease serine-like 1), KLK10, NES1 PRSSL1

**Immunogen:**

Recombinant Human Kallikrein-10 protein (31-274AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A