

# SAA1 Antibody, HRP conjugated



PACO33409

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Major acute phase reactant. Apolipoprotein of the HDL complex.

**Gene ID:**

SAA1

**Uniprot**

P0DJ18

**Synonyms:**

Serum amyloid A-1 protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA); Serum amyloid protein A(2-104); Serum amyloid protein A(3-104); Serum amyloid protein A(2-103); Serum amyloid protein A(2-102); Serum amyloid protein A(4-101)], SAA1

**Immunogen:**

Recombinant Human Serum amyloid A-1 protein (20-122AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A