

SAA1 Antibody, Biotin conjugated



PACO33411

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Major acute phase reactant. Apolipoprotein of the HDL complex.

Gene ID:

SAA1

Uniprot

P0DJ18

Synonyms:

Serum amyloid A-1 protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA); Serum amyloid protein A(2-104); Serum amyloid protein A(3-104); Serum amyloid protein A(2-103); Serum amyloid protein A(2-102); Serum amyloid protein A(4-101)], SAA1

Immunogen:

Recombinant Human Serum amyloid A-1 protein (20-122AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A