

# PLAU Antibody, Biotin conjugated



PACO33435

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.

**Gene ID:**

PLAU

**Uniprot**

P00749

**Synonyms:**

Urokinase-type plasminogen activator (U-plasminogen activator) (uPA) (EC 3.4.21.73)  
[Cleaved into: Urokinase-type plasminogen activator long chain A; Urokinase-type plasminogen activator short chain A; Urokinase-type plasminogen activator chain B],  
PLAU

**Immunogen:**

Recombinant Human Urokinase-type plasminogen activator protein (21-173AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A