

PACO33436

Product Information

Size:

50ug

Reactivity:

Helicobacter pylori

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the decarboxylation of orotidine 5'-monophosphate (OMP) to uridine 5'-monophosphate (UMP).

Gene ID:

pyrF

Uniprot

P56155

Synonyms:

Orotidine 5'-phosphate decarboxylase (EC 4.1.1.23) (OMP decarboxylase) (OMPDCase) (OMPdecase), pyrF

Immunogen:

Recombinant Helicobacter pylori Orotidine 5'-phosphate decarboxylase protein (1-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A