

ureA Antibody



PACO33440

Product Information

Size:

50ug

Reactivity:

Helicobacter pylori

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Ammonia produced by ureolysis increases the gastric pH thereby providing an environment permissive for colonization of the stomach.

Gene ID:

ureA

Uniprot

P14916

Synonyms:

Urease subunit alpha (EC 3.5.1.5) (Urea amidohydrolase subunit alpha), ureA, hpuA


Immunogen:

Recombinant Helicobacter pylori Urease subunit α protein (1-238AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ureA antibody at 2 μ g/ml + *Helicobacter pylori* bacteria liquid

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 27 kDa

Observed band size: 27 kDa

.