

PACO33476

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the biosynthesis of agmatine from arginine.

Gene ID:

speA

Uniprot

P21170

Synonyms:

Biosynthetic arginine decarboxylase (ADC) (EC 4.1.1.19), speA

Immunogen:

Recombinant Escherichia coli Biosynthetic arginine decarboxylase protein (3-656AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A