

# HLA-DPB1 Antibody



PACO33516

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000,  
IHC:1:20-1:200

**Protein Background:****Gene ID:**

HLA-DPB1

**Uniprot**

P04440

**Synonyms:**

HLA class II histocompatibility antigen, DP beta 1 chain (HLA class II histocompatibility antigen, DP(W4) beta chain) (MHC class II antigen DPB1), HLA-DPB1, HLA-DP1B

**Immunogen:**


Recombinant Human HLA class II histocompatibility antigen, DP & beta; 1 chain protein (30-225AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: HLA-DPB1 antibody at 2 $\mu$ g/ml

Lane 1: sw1990 whole cell lysate

Lane 2: HGC27 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO33516 at dilution of 1:50.