

HLA-DPB1 Antibody, Biotin conjugated



PACO33519

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HLA-DPB1

Uniprot

P04440

Synonyms:

HLA class II histocompatibility antigen, DP beta 1 chain (HLA class II histocompatibility antigen, DP(W4) beta chain) (MHC class II antigen DPB1), HLA-DPB1, HLA-DP1B

Immunogen:

Recombinant Human HLA class II histocompatibility antigen, DP β 1 chain protein (30-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A