

HLA-DPA1 Antibody, FITC conjugated



PACO33522

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HLA-DPA1

Uniprot

P20036

Synonyms:

HLA class II histocompatibility antigen, DP alpha 1 chain (DP(W3)) (DP(W4)) (HLA-SB alpha chain) (MHC class II DP3-alpha) (MHC class II DPA1), HLA-DPA1, HLA-DP1A HLASB

Immunogen:

Recombinant Human HLA class II histocompatibility antigen, DP α 1 chain protein (29-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A