

WNT3A Antibody, Biotin conjugated



PACO33619

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.

Gene ID:

WNT3A

Uniprot

P56704

Synonyms:

Protein Wnt-3a, WNT3A

Immunogen:

Recombinant Human Protein Wnt-3a protein (19-351AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A